

OVOL2 Antibody



PACO47206

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Zinc-finger transcription repressor factor. Plays a critical role to maintain the identity of epithelial lineages by suppressing epithelial-to mesenchymal transition mainly through the up-regulation of ZEB1 expression. Positively regulates neuronal differentiation. Suppresses cell cycling and terminal differentiation of keratinocytes by directly repressing MYC and NOTCH1.

Gene ID:

OVOL2

Uniprot

Q9BRP0

Synonyms:

Transcription factor Ovo-like 2 (hOvo2) (Zinc finger protein 339), OVOL2, ZNF339


Immunogen:

Recombinant Human Transcription factor Ovo-like 2 protein (19-114AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO47206 at dilution of 1:100.