

# KLHL3 Antibody, HRP conjugated



PACO47219

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Substrate-specific adapter of a BCR (BTB-CUL3-RBX1) E3 ubiquitin ligase complex that acts as a regulator of ion transport in the distal nephron. The BCR(KLHL3) complex acts by mediating ubiquitination of WNK4, an inhibitor of potassium channel KCNJ1, leading to WNK4 degradation.

**Gene ID:**

KLHL3

**Uniprot**

Q9UH77

**Synonyms:**

Kelch-like protein 3, KLHL3, KIAA1129

**Immunogen:**

Recombinant Human Kelch-like protein 3 protein (1-301AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A