

KLHL3 Antibody, Biotin conjugated



PACO47221

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Substrate-specific adapter of a BCR (BTB-CUL3-RBX1) E3 ubiquitin ligase complex that acts as a regulator of ion transport in the distal nephron. The BCR(KLHL3) complex acts by mediating ubiquitination of WNK4, an inhibitor of potassium channel KCNJ1, leading to WNK4 degradation.

Gene ID:

KLHL3

Uniprot

Q9UH77

Synonyms:

Kelch-like protein 3, KLHL3, KIAA1129

Immunogen:

Recombinant Human Kelch-like protein 3 protein (1-301AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A