

# CAMTA1 Antibody



PACO47258

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IF:1:50-1:200

**Protein Background:**

Transcriptional activator. May act as a tumor suppressor.

**Gene ID:**

CAMTA1

**Uniprot**

Q9Y6Y1

**Synonyms:**

Calmodulin-binding transcription activator 1, CAMTA1, KIAA0833

**Immunogen:**


Recombinant Human Calmodulin-binding transcription activator 1 protein (201-452AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunofluorescent analysis of HepG2 cells using PACO47258 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).