

ICE1 Antibody, FITC conjugated



PACO47264

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Component of the little elongation complex (LEC), a complex required to regulate small nuclear RNA (snRNA) gene transcription by RNA polymerase II and III. Specifically acts as a scaffold protein that promotes the LEC complex formation and recruitment and RNA polymerase II occupancy at snRNA genes in subnuclear bodies.

Gene ID:

ICE1

Uniprot

Q9Y2F5

Synonyms:

Little elongation complex subunit 1 (Interactor of little elongator complex ELL subunit 1), ICE1, KIAA0947

Immunogen:

Recombinant Human Little elongation complex subunit 1 protein (1213-1505AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A