

PACO47278

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytosol, adenosine deaminase activity, adenosine catabolic process, cellular response to DNA damage stimulus, hypoxanthine biosynthetic process, hypoxanthine salvage, inosine biosynthetic process, purine nucleotide interconversion, purine nucleotide salvage, purine ribonucleoside monophosphate biosynthetic process.

Gene ID:

add

Uniprot

P22333

Synonyms:

Adenosine deaminase (EC 3.5.4.4) (Adenosine aminohydrolase), add

Immunogen:

Recombinant Escherichia coli Adenosine deaminase protein (1-333AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A