

ADAM12 Antibody, HRP conjugated



PACO47283

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in skeletal muscle regeneration, specifically at the onset of cell fusion. Also involved in macrophage-derived giant cells (MGC) and osteoclast formation from mononuclear precursors.

Gene ID:

ADAM12

Uniprot

O43184

Synonyms:

Disintegrin and metalloproteinase domain-containing protein 12 (ADAM 12) (EC 3.4.24) (Meltrin-alpha), ADAM12, MLTN

Immunogen:

Recombinant Human Disintegrin and metalloproteinase domain-containing protein 12 protein (730-909AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A