

CD37 Antibody, Biotin conjugated



PACO47337

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

extracellular exosome, immunological synapse, integral component of plasma membrane, membrane, cell surface receptor signaling pathway.

Gene ID:

CD37

Uniprot

P11049

Synonyms:

Leukocyte antigen CD37 (Tetraspanin-26) (Tspan-26) (CD antigen CD37), CD37, TSPAN26

Immunogen:

Recombinant Human Leukocyte antigen CD37 protein (109-237AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A