

CTSE Antibody, FITC conjugated



PACO47364

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May have a role in immune function. Probably involved in the processing of antigenic peptides during MHC class II-mediated antigen presentation. May play a role in activation-induced lymphocyte depletion in the thymus, and in neuronal degeneration and glial cell activation in the brain.

Gene ID:

CTSE

Uniprot

P14091

Synonyms:

Cathepsin E (EC 3.4.23.34) [Cleaved into: Cathepsin E form I; Cathepsin E form II], CTSE

Immunogen:

Recombinant Human Cathepsin E protein (59-155AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A