

# EFEMP2 Antibody, Biotin conjugated



PACO47385

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

basement membrane, extracellular exosome, extracellular region, extracellular vesicle, extracellular matrix structural constituent.

**Gene ID:**

EFEMP2

**Uniprot**

O95967

**Synonyms:**

EGF-containing fibulin-like extracellular matrix protein 2 (Fibulin-4) (FIBL-4) (Protein UPH1), EFEMP2, FBLN4

**Immunogen:**

Recombinant Human EGF-containing fibulin-like extracellular matrix protein 2 protein (26-165AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A