

HNRNPLL Antibody, FITC conjugated



PACO47460

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

RNA-binding protein that functions as regulator of alternative splicing for multiple target mRNAs, including PTPRC/CD45 and STAT5A. Required for alternative splicing of PTPRC.

Gene ID:

HNRNPLL

Uniprot

Q8WVV9

Synonyms:

Heterogeneous nuclear ribonucleoprotein L-like (hnRNPLL) (Stromal RNA-regulating factor), HNRNPLL, HNRPLL SRRF

Immunogen:

Recombinant Human Heterogeneous nuclear ribonucleoprotein L-like protein (1-72AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A