

# HSPA9 Antibody



PACO47470

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Implicated in the control of cell proliferation and cellular aging. May also act as a chaperone.

**Gene ID:**

HSPA9

**Uniprot**

P38646

**Synonyms:**

Stress-70 protein, mitochondrial (75 kDa glucose-regulated protein) (GRP-75) (Heat shock 70 kDa protein 9) (Mortalin) (MOT) (Peptide-binding protein 74) (PBP74), HSPA9, GRP75 HSPA9B mt-HSP70

**Immunogen:**


Recombinant Human Stress-70 protein, mitochondrial protein (458-679AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO47470 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human breast cancer using PACO47470 at dilution of 1:100.