

HSPA9 Antibody, HRP conjugated



PACO47471

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Implicated in the control of cell proliferation and cellular aging. May also act as a chaperone.

Gene ID:

HSPA9

Uniprot

P38646

Synonyms:

Stress-70 protein, mitochondrial (75 kDa glucose-regulated protein) (GRP-75) (Heat shock 70 kDa protein 9) (Mortalin) (MOT) (Peptide-binding protein 74) (PBP74), HSPA9, GRP75 HSPA9B mt-HSP70

Immunogen:

Recombinant Human Stress-70 protein, mitochondrial protein (458-679AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A