

PLIN3 Antibody, Biotin conjugated



PACO47521

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Required for the transport of mannose 6-phosphate receptors (MPR) from endosomes to the trans-Golgi network.

Gene ID:

PLIN3

Uniprot

O60664

Synonyms:

Perilipin-3 (47 kDa mannose 6-phosphate receptor-binding protein) (47 kDa MPR-binding protein) (Cargo selection protein TIP47) (Mannose-6-phosphate receptor-binding protein 1) (Placental protein 17) (PP17), PLIN3, M6PRBP1 TIP47

Immunogen:

Recombinant Human Perilipin-3 protein (192-324AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A