

PEMT Antibody, Biotin conjugated



PACO47585

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes three sequential methylation reactions of phosphatidylethanolamine (PE) by AdoMet, thereby producing phosphatidylcholine (PC).

Gene ID:

PEMT

Uniprot

Q9UBM1

Synonyms:

Phosphatidylethanolamine N-methyltransferase (PEAMT) (PEMT) (EC 2.1.1.17) (EC 2.1.1.71) (PEMT2) (Phospholipid methyltransferase) (PLMT), PEMT, PEMPT PNMT

Immunogen:

Recombinant Human Phosphatidylethanolamine N-methyltransferase protein (2-49AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A