

RND3 Antibody, FITC conjugated



PACO47620

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Binds GTP but lacks intrinsic GTPase activity and is resistant to Rho-specific GTPase-activating proteins.

Gene ID:

RND3

Uniprot

P61587

Synonyms:

Rho-related GTP-binding protein RhoE (Protein MemB) (Rho family GTPase 3) (Rho-related GTP-binding protein Rho8) (Rnd3), RND3, ARHE RHO8 RHOE

Immunogen:

Recombinant Human Rho-related GTP-binding protein RhoE protein (152-239AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A