

# SCO1 Antibody, HRP conjugated



PACO47635

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Thought to play a role in cellular copper homeostasis, mitochondrial redox signaling or insertion of copper into the active site of COX.

**Gene ID:**

SCO1

**Uniprot**

O75880

**Synonyms:**

Protein SCO1 homolog, mitochondrial, SCO1, SCOD1

**Immunogen:**

Recombinant Human Protein SCO1 homolog, mitochondrial protein (29-176AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A