

# SLFN11 Antibody, Biotin conjugated



PACO47677

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Interferon (IFN)-induced antiviral protein which acts as an inhibitor of retrovirus protein synthesis. Specifically abrogates the production of retroviruses such as human immunodeficiency virus 1 (HIV-1) by acting as a specific inhibitor of the synthesis of retroviruses encoded proteins in a codon-usage-dependent manner. Binds to tRNAs and exploits the unique viral codon bias towards A/T nucleotides. The exact inhibition mechanism is unclear: may either sequesters tRNAs, prevents their maturation via post-transcriptional processing or accelerates their deacylation. Does not inhibit reverse transcription, integration or production and nuclear export of viral RNA. May play a role in cell cycle arrest and/or induction of apoptosis in response to exogenously induced DNA damage.

**Gene ID:**

SLFN11

**Uniprot**

Q7Z7L1

**Synonyms:**

Schlafen family member 11, SLFN11

**Immunogen:**

Recombinant Human Schlafen family member 11 protein (670-821AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A