

TLL1 Antibody



PACO47718

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Protease which processes procollagen C-propeptides, such as chordin, pro-biglycan and pro-lysyl oxidase. Required for the embryonic development. Predominant protease, which in the development, influences dorsal-ventral patterning and skeletogenesis.

Gene ID:

TLL1

Uniprot

O43897

Synonyms:

Tolloid-like protein 1 (EC 3.4.24), TLL1, TLL


Immunogen:

Recombinant Human Tolloid-like protein 1 protein (522-643AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO47718 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO47718 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human testis tissue using PACO47718 at dilution of 1:100.