

ZMPSTE24 Antibody, Biotin conjugated



PACO47757

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Proteolytically removes the C-terminal three residues of farnesylated proteins. Acts on lamin A/C.

Gene ID:

ZMPSTE24

Uniprot

O75844

Synonyms:

CAAX prenyl protease 1 homolog (EC 3.4.24.84) (Farnesylated proteins-converting enzyme 1) (FACE-1) (Prenyl protein-specific endoprotease 1) (Zinc metalloproteinase Ste24 homolog), ZMPSTE24, FACE1 STE24

Immunogen:

Recombinant Human CAAX prenyl protease 1 homolog protein (217-347AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A