

# LMNB1 Antibody, Biotin conjugated



PACO47797

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Lamins are components of the nuclear lamina, a fibrous layer on the nucleoplasmic side of the inner nuclear membrane, which is thought to provide a framework for the nuclear envelope and may also interact with chromatin.

**Gene ID:**

LMNB1

**Uniprot**

P20700

**Synonyms:**

Lamin-B1, LMNB1, LMN2 LMNB

**Immunogen:**

Recombinant Human Lamin-B1 protein (385-552AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A