

# PET100 Antibody



PACO47850

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200,  
IF:1:50-1:200

**Protein Background:**

**Gene ID:**

PET100

**Uniprot**

P0DJ07

**Synonyms:**

Protein PET100 homolog, mitochondrial, PET100, C19orf79

**Immunogen:**


Recombinant Human Protein PET100 homolog, mitochondrial protein (26-73AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human breast cancer using PACO47850 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO47850 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).