

# MEI4 Antibody



PACO47862

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:1000-1:5000

**Protein Background:**

Required for meiotic DNA double-strand break formation.

**Gene ID:**

MEI4

**Uniprot**

A8MW99

**Synonyms:**

Meiosis-specific protein MEI4, MEI4

**Immunogen:**


Recombinant Human Meiosis-specific protein MEI4-like protein (22-137AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Western Blot. Positive WB detected in: Mouse heart tissue. All lanes: MEI4 antibody at 2.7 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 45 kDa. Observed band size: 45 kDa.