

KLF5 Antibody, Biotin conjugated



PACO47877

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcription factor that binds to GC box promoter elements. Activates the transcription of these genes.

Gene ID:

KLF5

Uniprot

Q13887

Synonyms:

Krueppel-like factor 5 (Basic transcription element-binding protein 2) (BTE-binding protein 2) (Colon krueppel-like factor) (GC-box-binding protein 2) (Intestinal-enriched krueppel-like factor) (Transcription factor BTEB2), KLF5, BTEB2 CKLF IKLF

Immunogen:

Recombinant Human Krueppel-like factor 5 protein (69-346AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A