

TMEM240 Antibody



PACO47942

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

synaptic membrane.

Gene ID:

TMEM240

Uniprot

Q5SV17

Synonyms:

Transmembrane protein 240, TMEM240, C1orf70


Immunogen:

Recombinant Human Transmembrane protein 240 protein (42-87AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HepG2 cells using PACO47942 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).