

LINC00299 Antibody



PACO47970

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Gene ID:

LINC00299

Uniprot

Q6ZSB3

Synonyms:

Putative uncharacterized protein encoded by LINC00299, LINC00299, C2orf46
NCRNA00299


Immunogen:

Recombinant Human Putative uncharacterized protein encoded by LINC00299 protein
(2-105AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO47970 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).