

UNC13C Antibody



PACO48034

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

May play a role in vesicle maturation during exocytosis as a target of the diacylglycerol second messenger pathway. May be involved in the regulation of synaptic transmission at parallel fiber - Purkinje cell synapses.

Gene ID:

UNC13C

Uniprot

Q8NB66

Synonyms:

Protein unc-13 homolog C (Munc13-3), UNC13C


Immunogen:

Recombinant Human Protein unc-13 homolog C protein (414-645AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO48034 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO48034 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).