

# ATXN2 Antibody, HRP conjugated



PACO48083

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Involved in EGFR trafficking, acting as negative regulator of endocytic EGFR internalization at the plasma membrane.

**Gene ID:**

ATXN2

**Uniprot**

Q99700

**Synonyms:**

Ataxin-2 (Spinocerebellar ataxia type 2 protein) (Trinucleotide repeat-containing gene 13 protein), ATXN2, ATX2 SCA2 TNRC13

**Immunogen:**

Recombinant Human Ataxin-2 protein (174-383AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A