

PPP1R12C Antibody



PACO48114

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Regulates myosin phosphatase activity.

Gene ID:

PPP1R12C

Uniprot

Q9BZL4

Synonyms:

Protein phosphatase 1 regulatory subunit 12C (Protein phosphatase 1 myosin-binding subunit of 85 kDa) (Protein phosphatase 1 myosin-binding subunit p85), PPP1R12C, LENG3 MBS85


Immunogen:

Recombinant Human Protein phosphatase 1 regulatory subunit 12C protein (317-450AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human heart tissue using PACO48114 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO48114 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).