

CORO1C Antibody, Biotin conjugated



PACO48177

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in cytokinesis, motility, and signal transduction.

Gene ID:

CORO1C

Uniprot

Q9ULV4

Synonyms:

Coronin-1C (Coronin-3) (hCRNN4), CORO1C, CRN2 CRNN4

Immunogen:

Recombinant Human Coronin-1C protein (351-464AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A