

ABCG1 Antibody



PACO48206

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Transporter involved in macrophage lipid homeostasis. Is an active component of the macrophage lipid export complex. Could also be involved in intracellular lipid transport processes. The role in cellular lipid homeostasis may not be limited to macrophages.

Gene ID:

ABCG1

Uniprot

P45844

Synonyms:

ATP-binding cassette sub-family G member 1 (ATP-binding cassette transporter 8) (White protein homolog), ABCG1, ABC8 WHT1


Immunogen:

Recombinant Human ATP-binding cassette sub-family G member 1 protein (84-222AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO48206 at dilution of 1:100.