

ARHGDIB Antibody, Biotin conjugated



PACO48241

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Regulates the GDP/GTP exchange reaction of the Rho proteins by inhibiting the dissociation of GDP from them, and the subsequent binding of GTP to them. Regulates reorganization of the actin cytoskeleton mediated by Rho family members.

Gene ID:

ARHGDIB

Uniprot

P52566

Synonyms:

Rho GDP-dissociation inhibitor 2 (Rho GDI 2) (Ly-GDI) (Rho-GDI beta), ARHGDIB, GDIA2, GDID4, RAP1GN1

Immunogen:

Recombinant Human Rho GDP-dissociation inhibitor 2 protein (2-86AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A