

CHEK1 Antibody, Biotin conjugated



PACO48289

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CHEK1

Uniprot

O14757

Synonyms:

Serine/threonine-protein kinase Chk1 (EC 2.7.11.1) (CHK1 checkpoint homolog) (Cell cycle checkpoint kinase) (Checkpoint kinase-1), CHEK1, CHK1

Immunogen:

Recombinant Human Serine/threonine-protein kinase Chk1 protein (258-370AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A