

CNR2 Antibody, FITC conjugated



PACO48300

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Heterotrimeric G protein-coupled receptor for endocannabinoid 2-arachidonoylglycerol mediating inhibition of adenylate cyclase. May function in inflammatory response, nociceptive transmission and bone homeostasis.

Gene ID:

CNR2

Uniprot

P34972

Synonyms:

Cannabinoid receptor 2 (CB-2) (CB2) (hCB2) (CX5), CNR2, CB2A CB2B

Immunogen:

Recombinant Human Cannabinoid receptor 2 protein (305-360AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A