

CTRL Antibody, HRP conjugated



PACO48311

Product Information

Size:	Protein Background:
50ug	extracellular space, digestion, proteolysis.
Reactivity:	Gene ID:
Human	CTRL
Source:	Uniprot
Rabbit	P40313
Isotype:	Synonyms:
IgG	Chymotrypsin-like protease CTRL-1 (EC 3.4.21), CTRL, CTRL1
Applications:	Immunogen:
ELISA	Recombinant Human Chymotrypsin-like protease CTRL-1 protein (90-220AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A