

CTRL Antibody, Biotin conjugated



PACO48313

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

extracellular space, digestion, proteolysis.

Gene ID:

CTRL

Uniprot

P40313

Synonyms:

Chymotrypsin-like protease CTRL-1 (EC 3.4.21), CTRL, CTRL1

Immunogen:

Recombinant Human Chymotrypsin-like protease CTRL-1 protein (90-220AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A