

DRD1 Antibody, HRP conjugated



PACO48319

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Dopamine receptor whose activity is mediated by G proteins which activate adenylyl cyclase.

Gene ID:

DRD1

Uniprot

P21728

Synonyms:

D(1A) dopamine receptor (Dopamine D1 receptor), DRD1

Immunogen:

Recombinant Human D(1A) dopamine receptor protein (338-446AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A