

DUSP3 Antibody, FITC conjugated



PACO48324

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Shows activity both for tyrosine-protein phosphate and serine-protein phosphate, but displays a strong preference toward phosphotyrosines. Specifically dephosphorylates and inactivates ERK1 and ERK2.

Gene ID:

DUSP3

Uniprot

P51452

Synonyms:

Dual specificity protein phosphatase 3 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase VHR) (Vaccinia H1-related phosphatase) (VHR), DUSP3, VHR

Immunogen:

Recombinant Human Dual specificity protein phosphatase 3 protein (14-96AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A