

WNT7A Antibody, HRP conjugated



PACO48499

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. Signaling by Wnt-7a allows sexually dimorphic development of the mullerian ducts.

Gene ID:

WNT7A

Uniprot

O00755

Synonyms:

Protein Wnt-7a, WNT7A

Immunogen:

Recombinant Human Protein Wnt-7a protein (126-302AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A