

LEPR Antibody, HRP conjugated



PACO48515

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

LEPR

Uniprot

P48357

Synonyms:

Leptin receptor (LEP-R) (HuB219) (OB receptor) (OB-R) (CD antigen CD295), LEPR, DB OBR

Immunogen:

Recombinant Human Leptin receptor protein (874-1075AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A