

# LEPR Antibody, FITC conjugated



PACO48516

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

LEPR

**Uniprot**

P48357

**Synonyms:**

Leptin receptor (LEP-R) (HuB219) (OB receptor) (OB-R) (CD antigen CD295), LEPR, DB OBR

**Immunogen:**

Recombinant Human Leptin receptor protein (874-1075AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A