

FOXO1 Antibody, HRP conjugated



PACO48555

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FOXO1

Uniprot

Q12778

Synonyms:

Forkhead box protein O1 (Forkhead box protein O1A) (Forkhead in rhabdomyosarcoma), FOXO1, FKHR FOXO1A

Immunogen:

Recombinant Human Forkhead box protein O1 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A