

AQP6 Antibody



PACO48558

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Forms a water-specific channel that participates in distinct physiological functions such as glomerular filtration, tubular endocytosis and acid, base metabolism.

Gene ID:

AQP6

Uniprot

Q13520

Synonyms:

Aquaporin-6 (AQP-6) (Aquaporin-2-like) (Kidney-specific aquaporin) (hKID), AQP6, AQP2L


Immunogen:

Recombinant Human Aquaporin-6 protein (1-30AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO48558 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO48558 at dilution of 1:100.