

ODAPH Antibody



PACO48578

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Might promote nucleation of hydroxyapatite.

Gene ID:

ODAPH

Uniprot

Q17RF5

Synonyms:

Odontogenesis associated phosphoprotein, ODAPH, C4orf26


Immunogen:

Recombinant Human Odontogenesis associated phosphoprotein protein (24-130AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO48578 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human bladder cancer using PACO48578 at dilution of 1:100.