

# TGIF1 Antibody, FITC conjugated



PACO48612

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Binds to a retinoid X receptor (RXR) responsive element from the cellular retinol-binding protein II promoter (CRBP II-RXRE). Inhibits the 9-cis-retinoic acid, dependent RXR alpha transcription activation of the retinoic acid, responsive element. Active transcriptional corepressor of SMAD2. Links the nodal signaling pathway to the bifurcation of the forebrain and the establishment of ventral midline structures. May participate in the transmission of nuclear signals during development and in the adult, as illustrated by the down-modulation of the RXR alpha activities.

**Gene ID:**

TGIF1

**Uniprot**

Q15583

**Synonyms:**

Homeobox protein TGIF1 (5'-TG-3'-interacting factor 1), TGIF1, TGIF

**Immunogen:**

Recombinant Human Homeobox protein TGIF1 protein (8-375AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A