

# LIN28B Antibody, FITC conjugated



PACO48652

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

LIN28B

**Uniprot**

Q6ZN17

**Synonyms:**

Protein lin-28 homolog B (Lin-28B), LIN28B, CSDD2

**Immunogen:**

Recombinant Human Protein lin-28 homolog B protein (166-250AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A