

LIN28B Antibody, Biotin conjugated



PACO48653

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	LIN28B
Source:	Uniprot
Rabbit	Q6ZN17
Isotype:	Synonyms:
IgG	Protein lin-28 homolog B (Lin-28B), LIN28B, CSDD2
Applications:	Immunogen:
ELISA	Recombinant Human Protein lin-28 homolog B protein (166-250AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A