

PHIP Antibody, FITC conjugated



PACO48696

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable regulator of the insulin and insulin-like growth factor signaling pathways. Stimulates cell proliferation through regulation of cyclin transcription and has an anti-apoptotic activity through AKT1 phosphorylation and activation. Plays a role in the regulation of cell morphology and cytoskeletal organization.

Gene ID:

PHIP

Uniprot

Q8WWQ0

Synonyms:

PH-interacting protein (PHIP) (IRS-1 PH domain-binding protein) (WD repeat-containing protein 11), PHIP, WDR11

Immunogen:

Recombinant Human PH-interacting protein (508-608AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A