

# PHIP Antibody, Biotin conjugated



PACO48697

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Probable regulator of the insulin and insulin-like growth factor signaling pathways. Stimulates cell proliferation through regulation of cyclin transcription and has an anti-apoptotic activity through AKT1 phosphorylation and activation. Plays a role in the regulation of cell morphology and cytoskeletal organization.

**Gene ID:**

PHIP

**Uniprot**

Q8WWQ0

**Synonyms:**

PH-interacting protein (PHIP) (IRS-1 PH domain-binding protein) (WD repeat-containing protein 11), PHIP, WDR11

**Immunogen:**

Recombinant Human PH-interacting protein (508-608AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A