

HOXB13 Antibody, Biotin conjugated



PACO48701

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Gene ID:

HOXB13

Uniprot

Q92826

Synonyms:

Homeobox protein Hox-B13, HOXB13

Immunogen:

Recombinant Human Homeobox protein Hox-B13 protein (1-161AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A